

# HALYARD\* PURPLE NITRILE\* POWDER-FREE NON-STERILE AND STERILE NITRILE EXAM GLOVES

In any health care setting, protection from exposure to micro-organisms is of critical importance – for you and for the patients under your care. Gloves provide the primary shield required to prevent contact with infectious substances, including blood-borne pathogens, and reduce the risk of cross-contamination.

Our complete line includes gloves for primary care services, dental procedures, diagnostic tests and daily patient care.

The HALYARD\* PURPLE NITRILE\* product family offers a comprehensive line of latex-free non-sterile and sterile exam gloves. The PURPLE NITRILE\* is suitable for use across a wide variety of non-sterile medical applications. The PURPLE NITRILE\*-XTRA with a 5,5 cm longer cuff gives you extra protection in high risk situations involving exposure to blood borne pathogens or infectious fluids. The PURPLE NITRILE\* Sterile Pairs is excellent for sterile, non-surgical procedures in OB/GYN, ER/Trauma, ICU/CCU, wound care and burn centres.



## FUNCTIONAL BENEFITS

- Unique purple colour provides quick, visual differentiation from natural rubber latex and other synthetic gloves.
- Textured fingertips enhance wet and dry grip.
- Custom design enhances comfort and fit.
- Available in two lengths: 240 mm and 310 mm.
- Available in sterile pairs.

## CLINICAL SPECIFICATIONS

- Do not contain natural rubber latex: primary material is acrylonitrile-butadiene, reducing the potential for Type I glove-associated allergic reactions to natural rubber latex proteins.
- PURPLE NITRILE\*-XTRA material has been approved for use with a large amount of chemotherapy drugs and glutaraldehyde to determine permeation and breakthrough times.<sup>1</sup>
- Powder-free (according to ASTM D 6124), minimising the potential for powder-related complications, such as irritant contact dermatitis.

# PURPLE NITRILE\*

## Medical Gloves

### QUALITY STANDARDS

- Meets or exceeds AQL 1,5 in accordance with European Norm 455-1 and 455-2.
- CE-marked as Medical Device Class I (93/42/EEC).
- Biocompatible as tested according to ISO 10993 (1-2-5-10-12) which includes primary skin irritation, repeat challenge sensitisation and 200 person modified draize test.
- Sterilized in accordance with EN ISO 11137 Sterility Assurance Level (SAL) of  $10^{-6}$  (Sterile gloves).
- Manufactured in accordance with Quality System EN ISO 13485.

### PRODUCT SPECIFICATIONS

#### PURPLE NITRILE\*

MEASUREMENTS (MINIMUM)	GAUGE THICKNESS / MM (MEASURING SINGLE THICKNESS)	LENGTH
Middle Finger (15 mm from tip)	0,15	240 mm/310 mm (measured from the tip of the middle finger to the edge of the cuff)
Palm	0,12	
Cuff	0,09	

#### PHYSICAL PROPERTIES (minimum)

BEFORE AGING		AFTER AGING	
Force at Break	Ultimate Elongation	Force at Break	Ultimate Elongation
>9 N (Median 13-15 N)	> 500 %	>6 N (Median 13-15 N)	> 450 %

Note: The dimensions and properties listed above can vary within pre-established specifications. This document was created using the most recent information. In the interest of continuous improvement, the characteristics of the product may change without prior notice.

### ORDERING INFORMATION

#### PURPLE NITRILE\*

Length: 240 mm

PRODUCT	SIZE	GLOVES/BOX	BOXES/CASE
52000M	X-small	100	10
52001M	Small	100	10
52002M	Medium	100	10
52003M	Large	100	10
52004M	X-large	90	10

XTRA - Length: 310 mm

PRODUCT	SIZE	GLOVES/BOX	BOXES/CASE
50600M	X-small	50	10
50601M	Small	50	10
50602M	Medium	50	10
50603M	Large	50	10
50604M	X-large	50	10

Sterile Pairs - Length: 240 mm

PRODUCT	SIZE	PAIRS/BOX	PAIRS/CASE
52201M	Small	50	200
52202M	Medium	50	200
52203M	Large	50	200

References: 1. Test results on file. Warning: Do not use this product if you have a known allergy to chemical additives.



KNOWLEDGE NETWORK\* Clinical Education  
 Knowledgeable Customer Support  
 Expert Sales Force  
 Tools & Best Practices  
 Clinical Research  
 Commitment to Excellence

**For more information, please contact your Halyard Customer Care Team:**

Tel: 0800 917 65 85 (from UK & IE only)

Fax: 0800 169 02 35 (from UK & IE only)

customerservice.uk.ie@hyh.com

