

FLUIDSHIELD* Level 2

Surgical & Procedure Masks

ASTM F2100-11 Level 2 Compliant



COMFORTABLE FLUID RESISTANCE
AND SPLASH PROTECTION



FLUID-RESISTANT PROTECTION



For assured Level 2 fluid protection
look for this symbol on the package!

FLUIDSHIELD* Level 2

Surgical & Procedure Masks

ASTM F2100-11 Level 2 Compliant

CHOOSE THE MASK THAT FITS THE TASK

Designed to resist moderate fluid strike-through, new FLUIDSHIELD* Level 2 masks are available as either surgical masks with ties or procedure masks with earloops. If the fluid risk involves potential splash and splatter, both the surgical and the procedure mask styles are also available in a wraparound visor version for added protection.

Level 2 Surgical and Procedure Masks

Test Method	Test Description	ASTM F2100-11 Requirement Level 2	Level 2 Mask*
ASTM F1862	Fluid Resistance with synthetic blood	120 mm Hg	120 mm Hg
MIL-M-36954C	Delta P	< 5.0 mm H ₂ O	3.7 mm H ₂ O
ASTM F2101	Bacterial Filtration Efficiency (BFE)	≥ 98%	99.9%
ASTM F2299	Submicron Particle Filtration (PFE)	≥ 98% @ 0.1 micron	99.1%
16 CFR Part 1610	Flammability	Class 1	Class 1

*Based on the calculated median

EN14683

EU Standard Class	Bacterial Filtration Efficiency	Breathing Resistance (Pa)	Splash Resistance (kPa)	Microbial Cleanliness (cfu/g)
Type I	95 %	< 29.4	NA	≤ 30
Type II	98 %	< 29.4	NA	≤ 30
Type IIR	98 %	< 49.0	≥ 16.0	≤ 30

HALYARD* ASTM face masks are designed to meet industry standards and regulations, including EN14683. Some masks in the facial protection range designed for specific clinical applications meet EN 149:2001. FLUIDSHIELD* Level 2 masks comply with the requirements of Type IIR masks as defined by EN 14683.

LEVEL 2

AVAILABLE WITH OR WITHOUT WRAPAROUND VISORS.



Procedure Mask



Fog-Free Surgical Mask



Fog-Free Procedure Mask with visor



ASTM F2100-11
LEVEL 2



MODERATE FLUID PROTECTION

FLUIDSHIELD* Level 2 Surgical and Procedure Masks			
Item #	ASTM F2100-11 LEVEL 2 Product Description	Colour	Packaging
28810	Procedure Mask with earloops	Blue	50 masks/box 10 boxes/case
28821	Fog-Free Procedure Mask with wraparound visor, earloops	Blue	25 masks/box 4 boxes/case
28808	Fog-Free Surgical Mask with ties	Blue with cross hatch	50 masks/box 6 boxes/case
28809	Fog-Free Surgical Mask with wraparound visor, ties	Blue with cross hatch	25 masks/box 4 boxes/case

For more information, please contact your Halyard Customer Care Team:

UK & IE (toll-free)
Tel: 0800 085 82 04 - Fax: 0800 066 91 06
customerservice.uk.ie@hyh.com

Other countries
Tel: +32 2 585 51 14 - Fax: +32 2 897 00 95
customerservice.export@hyh.com

South Africa (toll-free)
Tel: 0801 111 332 - Fax: 0801 111 335
customerservice.za@hyh.com