

Cleanroom Nitrile Gloves

HALYARD* **PUREZERO*** NITRILE CLEANROOM GLOVES

HALYARD* offers an accelerator-free¹ portfolio of three non-sterile and two sterile cleanroom gloves to address the needs of pharmaceutical, medical device, microelectronics and semi-conductor manufacturing industries.

GLOVE SELECTION GUIDE

	Description	Designed for	Max Particle Count	ISO Class	Finish	Double Donning	Size Range	Fingertip Thickness	Case Count
NON-STERILE	HALYARD* PUREZERO* HG3 White Nitrile Gloves	semiconductor pharmaceutical medical device manufacturing	<950	Class 3 or higher	Tacky	Recommended for outer glove	XS-XL	.16mm (6 mil)	1000/cs
	HALYARD* PUREZERO* HG3 Light Blue Nitrile Gloves	semiconductor pharmaceutical medical device manufacturing	<1200	Class 3 or higher	Tacky	Recommended for outer glove	XS-XL	.10mm (4 mil)	1500/cs
	HALYARD* PUREZERO* HG3 SMOOTH White Nitrile Gloves	semiconductor pharmaceutical medical device manufacturing	<950	Class 3 or higher	Smooth	outer or under	XS-XL	.16mm (6 mil)	1000/cs
STERILE	HALYARD* PUREZERO* HG3 Light Blue Sterile Nitrile Gloves	pharmaceutical biotechnology sterile compounding aseptic processing	<1200	Class 3 or higher	Smooth	outer or under	6.0 6.5 7.0 7.5 8.0 8.5 9.0 10.0	.10mm (4 mil)	300 pairs/cs
	HALYARD* PUREZERO* HG3 White Sterile Nitrile Gloves	pharmaceutical biotechnology sterile compounding aseptic processing	<950	Class 3 or higher	Smooth	outer or under	6.0 6.5 7.0 7.5 8.0 8.5 9.0 10.0	.16mm (6 mil)	200 pairs/cs

Applies to all **PUREZERO*** Gloves:

Textured Fingertips

Accelerator-Free¹

AQL 1.5

For more information or samples, contact your distributor or visit: www.halyardhealth.com/cleanroom



¹ Not formulated with these commonly used vulcanizing chemicals: Sulfur, Thiurams, Thioxoles, Guanidines and Carbamates.
² Tested per ASTM D6319, EN 455-2
 *Registered Trademark or Trademark of O&M Halyard or its affiliates. ©2020. All rights reserved.
 H02757 H00078-12-01, 008271, COPY-03732